



Manufactured
in the UK



**Safelab Systems
Products**

Filtration Fume
Cupboards

Ducted Fume
Cupboards

Extract Systems

Cyclonic Scrubbers

Horizontal Laminar
Flow Cabinets

Vertical Laminar
Flow Cabinets

Safety Cabinets

Forensic Safety Cabinets

Polypropylene Cabinets

Glove Boxes

Storage Cabinets

**Safelab Systems
Service**

Maintenance Contracts

Airflow Testing

Ductwork Testing

Filter Testing

DOP Testing

KI Discus Testing

Installation

Commissioning

Lockable Service Docking Station for use with Airone 1000RS mobile fume cupboard



Standard docking station

New upgraded version
with RJ45 option.



**Designed specifically to offer an
additional safety feature within
the laboratory.**

**The (mobile fume cupboard)
service docking station stops unauthorised access to mains service
supplies.**

Water, waste and gas are securely locked away when not in use.

www.safelab.co.uk



Lockable Service Docking Station for use with Airone 1000RS mobile fume cupboard



Sturdy construction:
epoxy coated mild steel

Robust lock & strong hinges

Easily sited

Recess or flush mounted

Isolated 13amp single switch electrical socket

Quick release connections to water, waste and gas

Increased safety
within the laboratory.



Docking station with RJ45 option on
RHS of fume cupboard



Standard docking station on LHS of
fume cupboard

Safelab Systems Ltd

Unit 29 Lynx Crescent, Weston Super Mare, North Somerset. BS24 9BP

Tel: +44 (0) 1934 421 340

Fax: +44 (0) 1934 641 569

E-mail: safelab@safelab.co.uk

www.safelab.co.uk