



Recirculates clean air back into room



Airone BFSC-1000PP

Benchtop Forensic Science Safety Cabinet



The Airone BFSC-1000PP is a high performance Forensic Science Cabinet designed for the protection of evidence such as articles of clothing, mobile phones, wallets and bags from particulate contamination and dust.


It provides protection for the operator and the surrounding environment from noxious odours and other contaminants which may emanate from contents held within the cabinet.



Benefits

- Simple control panel
- Hour counter to display running time
- Visual indication/alarm of airflow
- Quiet operation/running mode
- Two speed fan
- Hanging rail
- Two Shelves
- White board: for evidence identification
- Lockable door
- Evidence lock
- Pre-filter

Technical Specification



- 1) Carbon and HEPA exhaust filters
- 2) Hanging rail
- 3) Mains Lead
- 4) Security tag location lugs
- 5) Door catches (1 lockable)
- 6) White board
- 7) Door operated fan speed switch
- 8) Illuminated green ON/OFF switch
- 9) Low airflow warning LED
- 10) Hour counter
- 11) Filter access panel/viewing window
- 12) Perforated shelves

Airone BFSC-1000PP dimensions:

External dimensions: (w x d x h) 1000 x 550 x 1300mm
Internal dimensions: (w x d x h) 975 x 520 x 880mm

(also available in 1200mm external width)

Option:

Internal lighting

Safelab manufacture a floor standing 1000mm or 1200mm wide Forensic Evidence Safety Cabinet with water/washdown facility



Contact:

Safelab Systems Ltd
Unit 29 Lynx Crescent
Weston Super Mare
North Somerset
BS24 9BP

Tel: 01934 421 340

Fax: 01934 641 569

E-mail: sales@safelab.co.uk

www.safelab.co.uk