



Manufactured
in the UK



**Safelab Systems
Products**

Filtration Fume
Cupboards

Ducted Fume
Cupboards

Extract Systems

Cyclonic Scrubbers

Horizontal Laminar
Flow Cabinets

Vertical Laminar
Flow Cabinets

Safety Cabinets

Forensic Safety Cabinets

Polypropylene Cabinets

Glove Boxes

Storage Cabinets

**Safelab Systems
Service**

Maintenance Contracts

Airflow Testing

Ductwork Testing

Filter Testing

DOP Testing

KI Discus Testing

Installation

Commissioning

Replacement filter
service

Airone BCSC-1000PP Benchtop Chemical Storage Cabinet



Designed specifically for the
safe storage of chemicals.

Protects staff and the
environment.

Benefits:

- Quiet operation/running mode
- Removable perforated shelves
- Large capacity carbon filter for the
removal of Formaldehyde

Optional internal lighting available

Dimensions:

External (WxDxH) 1000 x 550 x 1300mm
Internal (WxDxH) 980 x 525 x 880mm
Also available in 1200 width (external)

Construction & control system:

White polypropylene body & shelves
Clear PVC door

Mains On/Off

Two speed fan

Fan failure alarm

Safelab Systems Ltd

Unit 29 Lynx Crescent Weston Super Mare North Somerset BS24 9BP

Tel: +44 (0) 1934 421 340

Fax: +44 (0) 1934 641 569

E-mail: safelab@safelab.co.uk

www.safelab.co.uk