



Manufactured
in the UK



**Safelab Systems
Products**

Filtration Fume
Cupboards

Ducted Fume
Cupboards

Extract Systems

Cyclonic Scrubbers

Horizontal Laminar
Flow Cabinets

Vertical Laminar
Flow Cabinets

Safety Cabinets

Forensic Safety Cabinets

Polypropylene Cabinets

Glove Boxes

Storage Cabinets

**Safelab Systems
Service**

Maintenance Contracts

Airflow Testing

Ductwork Testing

Filter Testing

DOP Testing

KI Discus Testing

Installation

Commissioning

Replacement filter
service

Airone CSC-1000PP Chemical Storage Cabinet



Designed specifically for the safe storage of chemicals.

Protects staff and the environment.

Benefits:

Quiet operation/running mode

Removable perforated shelves

Large capacity carbon filter for the removal of Formaldehyde

Optional internal lighting available

Dimensions:

External (WxDxH) 1000 x 700 x 2140mm

Internal (WxDxH) 980 x 675 x 1530mm

Also available in 1200 width (external)

Construction & control system:

White polypropylene body & shelves

Clear PVC door

Mains On/Off

Fan failure alarm

Safelab Systems Ltd

Unit 29 Lynx Crescent Weston Super Mare North Somerset BS24 9BP

Tel: +44 (0) 1934 421 340

Fax: +44 (0) 1934 641 569

E-mail: safelab@safelab.co.uk

www.safelab.co.uk