



Manufactured
in the UK



**Safelab Systems
Products**

Filtration Fume
Cupboards

Ducted Fume
Cupboards

Extract Systems

Cyclonic Scrubbers

Horizontal Laminar
Flow Cabinets

Vertical Laminar
Flow Cabinets

Safety Cabinets

Forensic Safety Cabinets

Polypropylene Cabinets

Glove Boxes

Storage Cabinets

**Safelab Systems
Service**

Maintenance Contracts

Airflow Testing

Ductwork Testing

Filter Testing

DOP Testing

KI Discus Testing

Installation

Commissioning

Replacement filter
service

Airone CSC-1000PP Chemical Storage Cabinet



Designed specifically for the safe
storage of chemicals.

Protects staff and the
environment.

Benefits:

Quiet operation/running mode.

Removable perforated shelves.

Large capacity carbon filter for the
removal of a range of substance
dependent on your specific
requirements.

Optional internal lighting available.

Dimensions:

External (WxDxH) 1000 x 700 x 2140mm
Internal (WxDxH) 980 x 675 x 1530mm
Also available in 1200 width (external)

Construction & control system:

White polypropylene body & shelves
Clear PVC door

Mains On/Off

Fan failure alarm

Safelab Systems Ltd

Airone Building, 8 Beaufighter Road, Weston-Super-Mare BS24 8EE
Tel: +44 (0) 1934 421 340 Fax: +44 (0) 1934 641 569

E-mail: safelab@safelab.co.uk

www.safelab.co.uk