



Manufactured
in the UK



Safelab Systems Products

Filtration Fume
Cupboards

Ducted Fume
Cupboards

Extract Systems

Horizontal Laminar
Flow Cabinets

Vertical Laminar
Flow Cabinets

Safety Cabinets

Forensic Safety Cabinets

Polypropylene Cabinets

Storage Cabinets

Safelab Systems Service

Maintenance Contracts

Airflow Testing

Ductwork Testing

Filter Testing

DOP Testing

KI Discus Testing

LEV Testing

Installation

Commissioning

Lockable Service Docking Station

for use with mobile and
semi-mobile fume cupboards




**Allows quick and easy
connection of services utilising
bayonet connections.**

**The service docking station
stops unauthorised access to
mains service supplies, enabling
water, waste and gas to be securely
locked away when not in use.**




Mounting Options:

Can be mounted in three different ways. Options A and B using the standard docking station kit. C requires the optional docking station spacer (part number:



1. 13A switched socket outlet
2. 1/8" BSP natural gas coupling plug
3. 1/4" BSP water coupling plug
4. 1 1/2" BSP waste pipe coupler
5. Fume Cupboard tether attachment



Docking station on RHS of fume cupboard



Docking station on LHS of fume cupboard

Safelab Systems Ltd

Airone Building, 8 Beaufighter Road,
Weston Super Mare, North Somerset BS24 8EE

Tel: +44 (0) 1934 421 340

Fax: +44 (0) 1934 641 569

E-mail: safelab@safelab.co.uk

www.safelab.co.uk