



Manufactured  
in the UK



**Safelab Systems  
Products**

Filtration Fume  
Cupboards

Ducted Fume  
Cupboards

Extract Systems

Cyclonic Scrubbers

Horizontal Laminar  
Flow Cabinets

Vertical Laminar  
Flow Cabinets

Safety Cabinets

Forensic Safety Cabinets

Polypropylene Cabinets

Glove Boxes

Storage Cabinets

**Safelab Systems  
Service**

Maintenance Contracts

Airflow Testing

Ductwork Testing

Filter Testing

DOP Testing

KI Discus Testing

Installation

Commissioning

Replacement filter  
service

# Airone BCSC-1000PP Benchtop Chemical Storage Cabinet



Designed specifically for the  
safe storage of chemicals.

Protects staff and the  
environment.

**Benefits:**

Quiet operation/running mode.

Removable perforated shelves.

Large capacity carbon filter for the  
removal of a range of substances  
dependent on your specific  
requirements.

Optional internal lighting available.

**Dimensions:**

External (WxDxH) 1000 x 550 x 1300mm  
Internal (WxDxH) 980 x 525 x 880mm  
Also available in 1200 width (external)

**Construction & control system:**

White polypropylene body & shelves  
Clear PVC door

Mains On/Off

Two speed fan

Fan failure alarm

Safelab Systems Ltd

Airone Building, 8 Beaufighter Road, Weston-Super-Mare BS24 8EE

Tel: +44 (0) 1934 421 340

Fax: +44 (0) 1934 641 569

E-mail: safelab@safelab.co.uk

www.safelab.co.uk