



TECHNICAL SPECIFICATION

Manufacturer:	Safelab Systems Ltd
Address:	Airone Building, 8 Beaufighter Road, Weston-Super-Mare, BS24 8EE
Telephone:	+44 (0)1934 421340
E-mail:	enquiries@safelab.co.uk
Website:	www.safelab.co.uk
Supplier Contact:	Adrian Wyles – email: adrian@safelab.co.uk tel: +44 (0)7398 288292
Product description:	Lockable Service Docking station
Range:	Docking station for use with Airone RS mobile filtration & Airone DS Semi mobile ducted fume cupboards
Product Code:	Docking station
Application:	For use in Science Laboratories with Airone RS mobile filtration & Airone DS Semi mobile ducted fume cupboards. Designed for use in school laboratories. Allows quick and easy connection of services utilising bayonet connections. The docking station stops unauthorised access to mains service supplies, enabling water, waste, and gas to be securely locked away when not in use.
Standards:	Designed to comply with CLEAPSS G9 Fume cupboards in Schools.
Guarantee:	5 years warranty when Safelab carry out the annual service.
Installation:	Can be recess mounted or surface mounted – optional wall spacer is also available
External dimensions:	200mm wide x 110mm deep x 277mm high – excluding mounting bracket
External carcass:	Manufactured from 16swg mild steel and painted in epoxy powder coating
Lockable hinged access door:	Manufactured from 16swg mild steel and painted in epoxy powder coating – reversible – left or right-handed
Tether hook:	To accept gas safety stainless steel cable on the Airone RS & Airone DS
Weight:	4kg
Services:	
Gas fitting:	1/8" Quick release fitting - 1/8" BSPT to 10mm compression bulkhead fitting
Water fitting:	1/4" Quick release fitting – 1/4" BSPT to 10mm compression bulkhead fitting
Waste:	1 1/2" BSP waste pipe coupler & Blanking plug.
Electrical sockets:	Fascia mounted Single 13Amp outlet