





STAIRCASE AND STAIRWALKER TURNING CIRCLES AND CLEARANCE GUIDE



Clearance rotational requirement on staircases



Dimensions A and B are the minimum clearances required, allowing for any handrails (in blue).

You may be able to remove the wall-mounted handrails to create greater width.

FUME CUPBOARD TYPE	MINIMUM DIMENSION	
	A	B
Airone 1000RS	1100mm	1320mm
Airone 1000DS	1100mm	1320mm
Airone 1000X or 1000R	1100mm	1320mm
Airone 1200X or 1200R	1300mm	1475mm